INgene EHR

Prod Ref: 15.EHR.K.5

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I. TECHNICAL BASIS

In INgene EHR different techniques of molecular biology and immunology are combined in order to guarantee the reliability of the test and, as far as the complexity of this type of assay, allows its ease of performance:

1. Extraction techniques, specific retrotranscription and amplification techniques which guarantee the maximum degree of sensitivity and specificity.
2. Cloning techniques and design of plasmids which allow introducing exhaustive controls into each and every phase of the test. The INgene EHR incorporates external controls in a systematic way, which allows evaluating whether each of the assay steps has been carried out correctly. The kit includes an extraction, amplification and development control.
3. On-plate hybridisation techniques and techniques for labelling probes that allow the design of systems to read amplification results, which considerably increase the sensitivity of the test, facilitating it and – more importantly – “objectivising” the results interpretation process.

II. KIT COMPOSITION

			ELISA 48 reactions
		Volume	Nr Vials
AMPLIFICATION	Mixture A1 (specific primers for EHR)	750 µl	2
	Mixture B (enzyme mixture)	750 µl	2
	Positive control A1 – Amplification control	30 µl	1
DEVELOPMENT	Streptavidine Plate	48 poc.	1
	Development buffer	40 ml	1
	Desnaturalisation solution	15 ml	1
	Probe I – specific probe for EHR	90 µl	1
	Conjugate PO	90 µl	1
	Washing solution	250 ml	1
	Substrate	15 ml	1
	Stop solution	15 ml	1
	Positive control R1 (development) – ELISA Control	80 µl	1

To avoid any possible confusion, the reagents have been labelled with different colours depending on the process they are involved in. Thus, all reagents to be used during amplification have been labelled in **BLUE** and those used for developing in **YELLOW**.

III. GUIDELINES FOR THE CORRECT CONSERVATION OF THE KIT COMPONENTS

Supplied reagents must be maintained at +4°C.

IV. MATERIALS REQUIRED WHICH ARE NOT SUPPLIED

- VERY IMPORTANT: crushed ice to store the amplification reagents during handling.
- 0.2 or 0.5 ml Eppendorf tubes for thermocycler.
- Pipette tips with filter.
- 1 to 1000 µl pipettes.
- Assay tubes for the preparation of the amplification mixtures, probe solutions and conjugates for development, etc...
- Thermocycler.
- ELISA reader.
- Microfuge.

V. BASIC PRECAUTIONS TO AVOID CONTAMINATION

In the following, the precautions and recommendations to avoid any contamination and/or false results are detailed. If you are performing this assay for the first time, we recommend reading it through before using the kit. In general, the following precautions should be followed:

- Keep four separate areas for the preparation of samples, the preparation of the amplification mixtures, the amplification and the analysis of the products obtained (the latter two may overlap).
- In each designated area use material, lab-coats and equipment dedicated only to that area, i.e. do not allow it to cross boundaries.
- Change gloves every time any contamination is suspected.
- Never introduce amplified products into the sample preparation area or the area for the preparation of amplification mixture.
- Take special care when opening and closing the tubes. Avoid creation of aerosols.
- Maintain both reagents and samples closed as long as possible.
- Follow strict cleaning regulations inside the laboratory.

VI. COLLECTION AND TRANSPORT OF SAMPLES

The assay can be performed using any biological sample considered of interest by the clinician. Our technique has been standardised and contrasted with a considerable number of blood and serum samples. Other unrelated biological samples are not expected to limit the development of the technique.

The following are some recommendations and/or limitations relating to the method of obtaining and sending the samples to the laboratory:

- Use of heparin in any of the samples will limit the correct development of the assay, as it strongly inhibits the action of certain enzymes involved in the assay.
- Serum samples must be cooled from the moment they are obtained until they are processed. Processing needs to occur within 24 hours after obtaining the samples. In case a longer time is envisaged, we recommend freezing the samples.
- Tissue samples should be taken under the most aseptic conditions possible and be cooled immediately. As in the case of the sera, they will need to be kept refrigerated until they are processed, if the latter is scheduled within 24 hours, otherwise they should be frozen.

VII. EXTRACTION OF GENETIC MATERIAL

DNA EXTRACTION

The extraction of nucleic acids is a crucial point for the subsequent development of the PCR technique. Therefore, using a method that yields good-quality DNA is recommended (INGENASA Extraction could be provided under request).

Determining the quantity and quality of the DNA

Once the DNA has been isolated, it is recommended to check its quality and quantity.

SUGGESTED PROTOCOL

REQUIRED MATERIALS (Not supplied)

- Extraction buffer.
- DNAase-RNAase-free water.
- Chloroform.
- Isopropanol.
- Ethanol 70% (-20°C).
- Sterile PBS (only in the case of tissue samples).
- 1.5 ml Eppendorf tubes.
- Microfuge.

PROCEDURE

Initial preparation of tissue samples:

Cut a piece of tissue of approx. 2 mg and add sterile PBS until a concentration of approximately 0.4 mg/ml is reached. Homogenise the mixture under sterile conditions using a syringe or any other method ensuring complete homogenisation.

The samples treated in this way can follow the same extraction procedure as the serum samples.

1. Prepare 1.5 ml Eppendorf tubes for each of the samples to be analyzed plus one additional tube for negative extraction controls.
2. Add 100 µl of sample (serum or processed tissue according to previous procedure) and 750 µl extraction buffer. Follow the same procedure for the extraction control.
3. Add 0.25 ml of chloroform per every 0.75 ml of extraction buffer. Shake the tubes manually for 15 seconds.
4. Homogenise the mixtures and keep them at room temperature for 5 minutes.
5. Centrifuge the tubes at 12000 g for 15 minutes at +4°C.
6. Two phases can be differentiated after centrifugation, an aqueous one in the upper part and another one with an orange tinge, in the bottom

- part of the tube. The DNA remains in the upper aqueous phase.
- Prepare as many 1.5 ml Eppendorf tubes as samples being processed, label the tubes and add 0.5 ml of isopropanol to each one. Transfer the aqueous phase, approx. 350 µl, to the tubes previously prepared with isopropanol. Homogenise each mixture and incubate for 10 minutes at -20°C.
 - Centrifuge again for 10 minutes at 12000 g and decant the supernatant. The DNA stays in the pellet.
 - Add 1 ml of ethanol 70% (cold, stored at -20°C) to the tube and mix carefully.
 - Centrifuge at 7500 g for 15 minutes at +4°C.
 - Remove the supernatant by letting the excess liquid drain thoroughly on top of a filter paper. Leave the tubes open for 10 minutes to allow them to dry. **CAUTION:** *Never allow them to dry completely as the DNA will become insoluble if the pellet dries out completely.*
 - Dissolve the DNA in 50 µl of DNase-free water, thoroughly homogenise the mixture and store at 4°C until processing. **CAUTION:** *if amplification is not immediate, keep DNA at -20°C.*

VIII. AMPLIFICATION OF THE GENETIC MATERIAL

REQUIRED MATERIALS

- Crushed ice.
- DNA extracted from the samples.
- Mixture A1 – **KEEP IN CRUSHED ICE AT ALL TIMES**
- Mixture B – **KEEP IN CRUSHED ICE AT ALL TIMES WHILE OUTSIDE THE FREEZER.**
- Amplification control – Control A1.
- DNase-RNase-free water.

PROCEDURE

- Prepare and appropriately label as many tubes for the amplification as samples to be processed adding an additional tube for the negative extraction control.
 - Take mixture A1 and mixture B out of the freezer keeping them in crushed ice. Shake them lightly to ensure correct homogenisation. Prepare an appropriate amount of amplification mixture for the number of samples to be processed. The volume of each reagent to be mixed for each of the samples is:

Mixture A1	25 µl
Mixture B	25 µl
Final volume of mixture	50 µl
 - The tube used for mixing should be kept in crushed ice at all times. Likewise, it is recommendable to prepare an excess amount of mixture (calculate an extra 10% for all reagents) in order to compensate for possible losses of volume during pipetting. Once the mixture is prepared, homogenise correctly. Place the tubes previously labelled in crushed ice and add 50 µl of the mixture prepared in this way to each tube.
 - Add 3 µl of previously extracted DNA samples to each tube, 3 µl of positive control A1 (amplification control) to the corresponding tube and 3 µl of water to the tube labelled as negative control. Carefully mix the contents of each tube and ensure that all the liquid has been deposited at the bottom of the tube. If not, the tubes may be centrifuged lightly until this occurs.
 - Set the thermocycler to the following conditions:
 - **Pre-denaturalisation:** 3 min at 94 °C
 - **PCR Amplification:**

1 min at 94°C	}	40 Cycles
1 min at 72°C		
 - **Final extension:** 1 cycle of 5 min. at 72°C
- Maintain samples at 4°C until removal from thermocycler.

IX. DETECTION ELISA

MATERIALS

- Streptavidin plates
- Development Positive Control (R1)
- Developing buffer
- Denaturalisation solution – **CAUTION:** This reagent is a strong base and must be handled with care. This reagent is four times concentrated and needs to be diluted at the time of use.
- Probe I – specific for EHR
- Conjugate peroxidase (PO) – 100 x concentrated. Dilute conjugate only immediately before use and only prepare exact volume required, as it remains active for a short time once diluted.
- Washing solution.
- Substrate
- Stop solution

ASSAY PERFORMANCE PROTOCOL

Technical basis of the ELISA on-plate development.

The oligonucleotides used in the INGene EHR assay have been previously labelled with biotin so that the

amplification products – in case amplification occurs – will link to the streptavidin which lines the wells of the plate.

Following denaturalisation of the linked fragments, we add a specially designed probe which is complementary to an interDNA region of the amplified fragment that has also previously been labelled.

Under favourable conditions, hybridisation will occur in those wells where positive samples are being processed, and this hybridisation may be revealed by adding a PO conjugate which is able to specifically link the labelled probe. This in turn, can easily be developed by adding the adequate substrate to the PO as in any ELISA.

Procedure

- Preparation of the plate

Take out the required number of wells from the fridge and allow standing for at least 15 minutes before starting the assay.

You will need:

- 2 wells for the ELISA controls (one for the positive and one for the negative).
- 1well for each sample analyzed which – apart from the actual sample – includes all extraction and amplification controls which have been introduced into the assay.

Conveniently label each of the wells.

2.- Adding the samples

Once they reach room temperature, dispense 150 µl of developing buffer into each of the wells. Add 10µl of the sample to be analyzed (the product of the PCR) to each well. Homogenise the content of each well taking care to avoid any cross-contamination of wells, and then incubate the plate for 30 minutes at room temperature.

3.- Denaturalisation of the DNA chains

Prepare 1x denaturalisation solution from the one supplied (4x) in the following way:

$$Vf = n \times 170\mu l$$

Denaturalisation solution 1x = Vf/4 of denaturalisation solution + 3Vf/4 of H₂O.

Without removing the content, add 150 µl of denaturalisation solution to each well. Carefully homogenise the mixture and maintain at room temperature for 10 minutes (if possible, shake during incubation).

4.- Washing

Wash the wells 5 times with washing solution.

RECOMMENDATIONS FOR WASHING CYCLES:

The washing cycles can be performed using an automatic plate washer or a micropipette which allows dispensing 300 µl per well. Following incubation, wash according to the following instructions:

Remove the content of the plate by swiftly turning it upside down so as to avoid the exchange of fluids between wells. As a precaution, emptying of the wells should be done with care, as the denaturalising solution is highly irritating.

Dispense 300 µl of washing solution per well.

Carefully shake the plate avoiding the exchange of material between wells.

Swiftly turn the plate upside down to empty its content.

Repeat the process as often as indicated in the procedure.

Before removing the content of the last wash, make sure the next reagent is ready for immediate use. The plate must not stay dry.

After the last wash shake the plate upside down over an absorbent filter paper.

5.- Adding the specific probe

They are supplied 100x concentrated; therefore they need to be diluted in developing buffer before use. Volume of diluted probe to be prepared should be the exact volume required for the number of wells to assay, as any leftover will have to be discarded. The volume required in each case is:

$$(1.7 \mu l \text{ of probe in } 170 \mu l \text{ of developing buffer}) \times n$$

Where n is the number of wells to which the specific probe will be added.

Incubate the plate for 30 minutes at 57°C.

Remove the content from the wells.

6.- Wash

Wash the wells 5 times following the procedure specified under step 4.

7.- Adding of the conjugate

The conjugate is supplied 100x concentrated; therefore it needs to be diluted accordingly in developing buffer before use. As this reagent easily loses its activity once diluted, prepare only the exact volume required for the number of wells being used. The required volume is as follows:

$$(1.7 \mu l \text{ of conjugate supplied in } 170 \mu l \text{ of developing buffer}) \times n$$

Where n is the total number of wells in use

Once prepared, add 150µl to each of the wells and incubate for 20 minutes at room temperature.

8.- Wash

Wash the wells 4 times with washing solution according to the procedure specified under step 4, and wash once more with distilled water. Drain the plate carefully over a filter paper.

9.- Development

Add 150 µl of substrate to each well and allow the plate to stand at room temperature and in darkness for 10 minutes. Then add 100µl of stop solution.

10.- Reading

Read the OD values for each well with a spectrophotometer using a filter of 450 nm.

Validation of the assay

Before proceeding to the interpretation of the results obtained in the amplification assay, one has to verify that the ELISA has been performed correctly. For this, the ELISA control values need to be checked:

- The OD value for the negative ELISA control must be lower than 0.220
- The OD value for the positive ELISA control must be higher than 1.

If these values are not fulfilled, the ELISA procedure has not been performed correctly and needs to be repeated. THIS DOES NOT IMPLY THAT THE PCR HAS BEEN CARRIED OUT INCORRECTLY.

Evaluation of the samples

Samples which develop OD values above 0.155 units of the OD value of the negative control are considered **POSITIVE**.

Samples which develop OD values below this value are considered **NEGATIVE**.

X. INTERPRETATION OF RESULTS

Validation of the extraction procedure

To check whether the extraction process has been completed correctly, the behaviour of the extraction controls in the ELISA is evaluated.

If the positive extraction control can be considered positive in the ELISA while the negative extraction control stays within the range of negative values, the extraction process can be regarded as successfully completed.

control stays within the range of negative values, the amplification process can be regarded as successfully completed.

Interpretation of results

Once the possibility of errors in the process has been excluded, the result of the samples can be evaluated.

Validation of the amplification procedure

To check whether the amplification process has been completed correctly, the behaviour of the amplification controls in the ELISA is evaluated.

If the positive amplification control can be considered positive in the ELISA while the negative amplification

All the samples which show positive OD results in the ELISA should be considered positive, i.e. they contain genetic material of EHR.

All the samples that have shown negative OD values in the ELISA should be considered **NEGATIVE**, i.e. they do not contain genetic material of EHR.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.es
www.ingenasa.es