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
POSITIVO

I. TECHNICAL BASIS

Ingezim Heart-crom is a rapid qualitative one step immunoassay using sandwich method with monoclonal antibodies directed against circulating *Dirofilaria immitis* antigen.

II. STORAGE:

All the reagents must be kept between 4 and 25 °C.

IMPORTANT! DO NOT FREEZE

III. PRECAUTIONS:

1. Read carefully the instructions of use.
2. **Kept the reagents at room temperature before their use.**
3. Do not mix reagents or instruction from different kits.
4. Avoid any contamination of the reagents.
5. Do not use the kits after the expiry date.
6. There should be no eating, drinking or smoking where specimens or kit reagents are being handled.
7. Do not pipette by mouth.

IV. TEST PROCEDURE:

1. Preparing the test device:

Removed the test device from the protective pouch and place it on a flat surface.

NOTE: Do not remove the device from the pouch until ready to use.

2. Sample addition:

Gently mix the sample by inverting. Using the provided pipette, dispense two drops of sample (serum, plasma or whole blood) into window n° 1. Wait for the sample to be absorbed.

Serum or plasma: Fresh, previously frozen or stored at 2°-7° C, may be used in this test. Serum or plasma may

be stored for up to 7 days at 2°-7° C. For longer storage, samples should be frozen.

Whole blood: Must be anticoagulated (EDTA, heparin, etc.) and may be used either fresh or after refrigeration for up to 5 days. Hemolyzed sample will not affect results. Previously frozen or older samples must be centrifuged before use.

PRECAUTION: *Never touch the nitro-cellulose membrane with the pipette.*

3. Solution buffer addition:

Holding the diluent bottle vertically, promptly add 4 drops of the buffer to the window n°1.

4. Read and interpretation:

After 5 minutes read the results; do not read the results later than 15 minutes:


POSITIVE RESULTS: the test is positive if two coloured bands appear. One band will appear in the window n° 2 (Test line) and one in the window n°3 (Control line). Any coloured band in the window n° 2 (test line) should be considered positive.

NEGATIVE RESULT: There is no line in the well "2" (result well) and there is a purple coloured line in the well "3" (control well).

NO VALID TEST: If there is no line into the well "3" the assay can't be considered valid. If this occurs, add one additional drop of diluent and wait for 5 minutes. If a coloured band doesn't appear in the window n°3, the test is invalid and should be repeated.

Positive results may be observed in as short as 1 minute depending on the severity of the infection. **Coloured bands, which appears after 15 minutes are not diagnostic, and should be ignored.**

Add 2 drops of sample and 4 drops of buffer solution into the well "1"

Wait 5 minutes and read:

NEGATIVE

POSITIVE

Developed and manufactured in
Spain by:

Distributed in Spain by:

INMUNOLOGIA Y GENETICA APLICADA,
S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

IT-73840 ISO 14001:2004 ISO 9001:2008
IT-73780 9191.INGE 9175.ING2