INGEZIM LEISHMACROM (Anticuerpo)

Prod Ref: 15.LSH.K41

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

POSITIVO

I. TECHNICAL BASIS

The technical basis of Ingezim Leishmacrom is an indirect immunochromatography, where we are going to detect the presence of specific antibodies to Leishmania.

II. STORAGE:

All the reagents must be kept between 4 and 25 °C.

IMPORTANT! DO NOT FREEZE

III. PRECAUTIONS:

1. Read carefully the instructions of use.
2. **Kept the reagents at room temperature before their use.**
3. Do not mix reagents or instruction from different kits.
4. Avoid any contamination of the reagents.
5. Do not use the kits after the expiry date.
6. There should be no eating, drinking or smoking where specimens or kit reagents are being handled.
7. Do not pipette by month.

IV. TEST PROCEDURE:

1. Preparation of the sample:

The sera sample needs to be at 1/50 dilution, so take 10 µl of the sera sample and put it into a diluent tube.

Shake for getting a right mixture.

2. Sample addition:

Dispense 10 µl of the sample dilution previously made in the lower section of the window "3" NOTE: Open the aluminium bag just in the moment to make the assay.

PRECAUTION: *The pipette never must touch the nitro-cellulose membrane.*

Chromatographic Buffer addition:

1 minute later add the chromatographic buffer solution (5 drops) in the window "1".

3. Read the results:

Read the results after 10 minutes.


- 1 In order to be sure that the chromatography has run properly, a

pink line must appear at any case in the window nº 3 (control line window)


- 1 If on the window nº 2 (test window), appears a pink line the sample assayed must be reported as **POSITIVE** to Leishmania antibodies. It doesn't matter the intensity on the colour, a presence of line in this window means presence of antibody in the sample
- 1 If there is no line in the window "2" (result window) and there is a pink coloured line in the window "3" (control window) the sample must be reported as **NEGATIVE** to Leishmania antibodies.
- 1 **NO VALID TEST:** If there is no line into the window "3" the assay can't be considered valid.

No result must be considered before 20 minutes. Any line arising after this period of time must be not considered.


1. - Preparation of the sample:

2. - Add 10 μ l of the sample dilution **on the lower part of window "3"**.

3. - After 1 minute add 5 drops of the buffer solution in the window "1".

4. - Wait **10 minutes** for reading results

NEGATIVE

POSITIVE

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by:

IT-73840 ISO 14001:2004 ISO 9001:2008
IT-73780 9191.INGE 9175.ING2