## INGEZIM FIV-VET

Prod Ref: 16.FIV.K8-32

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## I. TECHNICAL BASIS

---

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on plate. After washing to eliminate all non-fixed material from the sera sample, we can detect the presence of cat immunoglobulins using a specific peroxidase conjugate.

After addition of the substrate a colorimetric reaction will appear which could be detected

In this way the presence of colour means the presence of antibodies against the virus in the cat sera, and the absence of colour the absence of specific antibodies.

The aim of this kit is to provide the users with a reliable diagnostic technique for this disease.

In our kit is important to remark the use of a synthetic peptide like antigen. This method warrants the total absence of infectivity in the kit.

## II. PRECAUTIONS AND WARNING FOR USERS:

---

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.the pipette.
8. Use a new loop for each sera sample.

## III. STORAGE CONDITIONS:

---

All reagents and components, must be stored at 4 °C.

## IV. INFORMATION ABOUT THE WASHING STEPS:

---

The washing steps could be done using a washing solution bottle or a multichanel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instruction:

1. Throw out the content of the wells by a brusque turn over of the to avoid the possible mixture of the content from one well to another
2. Fill up the wells with the washing solution (prepared following the instructions).
3. Shake delicately the wells, avoiding the contamination between wells.
4. Turn over the plate brusquely to empty the wells.
5. Repeat the process as many times is indicated on the instructions of the kit.
6. Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
7. Do not maintain the plate on dry more time than strictly needed.
8. After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

The sera samples need to be at **1/100** dilution in sera diluent to be tested.  
You get it using the calibrated loop and following the instruction

## VI. WAY TO USE THE CALIBRATED LOOPS:

- Ⓡ Use one different calibrated loop for each sera sample.
- Ⓡ For taking the adequate volume introduce only the final loop into the sample container.
- Ⓡ The captured volume (1 µl) is deposited by simple introduction of the loop into the appropriate well (in which previously it has been added 2 drops of sera diluent).
- Ⓡ The dilution is homogenised by soft turns of the loop into the well.

## VII. PREPARATION OF REAGENTS:

All the reagents are ready to use, without any previous manipulation. Only the washing solution is 10 x concentrated.

distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

### WASHING SOLUTION PREPARATION:

Dilute one part of the concentrate washing solution provided in the kit into 9 parts of

**Example:** 60 ml. of concentrated washing solution into 540 ml. of distilled or deionized water.

## VIII. TEST PROCEDURE:

1. Take out the strips or wells to be used and keep it at room temperature before starting the test (at least 2 hours).
2. Add **2 drops** of serum diluent on every well where serum samples are going to be tested. Then and using the calibrated loops add **1 µl** of serum sample on each well following the previous instructions. Positive and negative controls are supplied already diluted, that's why it is only needed to add **2 drops** of each sera control into two different empty wells. Seal the wells and incubate for **10 minutes at room temperature**.
3. **Wash 5 or more times** following the described procedure.
4. Add **2 drops** of conjugate, to each well. Seal the plate and incubate for **10 minutes at room temperature**.
5. **Wash 5 or more times** following the described procedure.
6. Add **1 drop** of chromogen solution, to each well and **1 drop** of substrate solution in the same way. Keep the plate for **5 min at room temperature**.  
Read the results.

## IX. READING AND INTERPRETAION OF THE RESULTS:

The reading must be done by visual colour turn detection.

**Positive control:** There must be a colour development (blue).

### Kit validation:

The colour development in the wells of the control sera, will be used as reference:

### Results Interpretation:

**Negative samples:** There is no colour turn detection.

**Negative control:** There must be no colour development (transparent).

**Positive samples:** There is colour turn detection.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in

by:

IT-73840  
IT-73780

ISO 14001:2015  
9191.INGE

ISO 9001:2015  
9175.ING2