INGEZIM ASFV-CSFV CROM Ab

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
I. TECHNICAL BASIS

This immunochromatographic test uses proteins of African swine fever virus (ASFV) and Classical swine fever virus (CSFV) to detect specific antibodies in blood or serum samples.

The test is based on the use of three microspheres: blue latex beads which are covalently linked to E2 protein of CSFV, red latex beads which are linked to VP72 protein of ASFV, and green latex beads which are linked to a control protein.

On the membrane, two test lines are printed: T1 formed by the E2 protein of CSFV and T2 formed by the VP72 protein of ASFV. Furthermore, the test contains a control line formed by a specific monoclonal antibody (MAb) for the control protein.

The antibodies present in the sample react with the latex particles which are coated with specific target proteins. This *latex particles/protein/antibodies* complex flows through the membrane. Depending on the type of antibody present in the sample, different coloured lines appear, which will be used to interpret the test results.

II. STORAGE

All the reagents must be kept between +4°C and +25 °C.

1. PRECAUTIONS *Read the use instructions carefully.*

2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration date and do not mix components from different batches.
6. Do not eat, drink, or smoke where specimens or kit reagents are being handle.

III. SAMPLES

The test has been designed for use with serum or blood samples. Serum samples must be fresh, refrigerated up to 7 days at 2-8 °C or stored at - 20 ° C. Blood samples must be fresh or refrigerated up to 4 days at 2-8 °C, and collected with any anticoagulant (EDTA, heparin, etc.). Blood samples without any anticoagulant, could contain micro-clots that can block the device and result in nonspecific reactions.

IV. TEST PROCEDURE

- **Allow all kit components and samples to reach room temperature, before use.**
- Take 20 µl of the serum sample, add it in a clean well.
- Add 3-4 drops (100 µl) of running buffer on the same well.
- Take the strips out of the tube and put them vertically, one in each well containing the prepared samples as describe previously. The arrows should be pointing at the sample direction.
- Keep the strips in the wells during 10 minutes. After these 10 minutes, take the strips out and put them on a flat and (preferably) white surface, for reading the results.

IT IS VERY IMPORTANT TO FOLLOW THESE INSTRUCTIONS AND RESPECT THE INDICATED VOLUMES TO OBTAIN VALID RESULTS

V. INTERPRETATION OF RESULTS



Interpretation of results should be done as follows, 10 minutes after the addition of the sample:

NEGATIVE: A single GREEN line appears in the result area. This line must always appear.

POSITIVE: ASFV: A GREEN line and a second PINK/RED line appear in the result area. Colour intensity may vary depending on the antibody concentration. CSFV: A GREEN line and a second BLUE line appear in the result area. Colour intensity may vary depending on the antibody concentration.

INVALID: If no green line appears, the test will be invalid. Repeat the test using a new device.

- **No results must be considered after 10 minutes.**

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENÉTICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91 368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: