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## I. TECHNICAL BASIS

---

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a solid support (polystyrene plate). When a serum sample contains specific antibodies against the virus, they will bind the antigen adsorbed on plate. After washing to eliminate all non-

fixed material from the sera sample, we can detect the presence of horse immunoglobulins using a specific peroxidase conjugate. After the addition of the substrate a colorimetric reaction will appear which could be measured by an ELISA reader.

In this way the presence of colour means the presence of antibodies against the virus in the horse sera, and the absence of colour the absence of specific antibodies.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Substrate must be handled with care, it is very sensible to light and contamination.
11. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water

## III. STORAGE OF COMPONENTS

---

All reagents and plates must be stored between +2°C and +8°C.

The stop solution is precipitated at this temperature. Incubate at 37°C before using, to eliminate the precipitation.

## IV. INFORMATION ABOUT THE WASHING STEPS

---

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

---

Serum samples must be diluted 1/100 in provided diluent (i.e, 5 µl of sample and 495 µl of diluent).

## VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled water. When ready, this solution remains stable between +2°C y +8°C.

- **Diluents:**

Dilute 1 volume of provided 5x concentrated diluent with 4 volumes of distilled or deionized water.

**Controls:**

Controls are ready to use.

- **Preparation of the conjugate:**

Prepare the dilution 1/100 in provided diluent (for one plate, 110 µl of concentrated conjugate up to 11 ml of diluent)

**ATTENTION:** For conjugate, prepare only the quantity needed for each time because the remainder volume has to be rejected.

## VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of each sample per well of the plate (diluted as previously indicated on point V). We recommend duplicating each sample for confirmatory purposes (not mandatory).  
  
Add 100 µl of the positive and negative controls, to their respective wells (two wells/control). Seal the plate and **incubate for 60 min. at 37°C**.
3. Wash 4 times following the described procedure.
4. Add 100 µl of conjugate (prepared following previous instructions) to each well. Seal the plate and **incubate for 30 min. at 25°C**.
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate for **15 min at room temperature**. We recommend using a multichannel pipette to make the process quickly.
7. Add 100 µl of stop solution to each well to stop the reaction. We recommend adding this reagent following the same order as it was added the substrate.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm**.

### A. Validation criteria:

The test is considered valid when:

- OD value of the positive control is higher than 1.0
- OD value of the negative control is lower than 0.3

### B. Positivity index calculation:

If you run samples in duplicate, OD values will be calculated for each sample as the arithmetic mean of both values.

To obtain S/P of the samples:

$$S/P = \frac{OD\ 405\ nm\ Sample - OD\ 405nm\ C\ (-)}{OD\ 405nm\ C\ (+) - OD\ 405nm\ C\ (-)}$$

### C. Results Interpretation (Cut off):

When you are running duplicate samples, OD values will be calculated for each sample as the arithmetic mean of both values

- Samples with S/P ≥ 0.3 will be considered as positive.
- Samples with S/P < 0.3 will be considered as negative.

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in \_\_\_\_\_ by:

IT-73840  
IT-73780

ISO 14001:2015  
9191.INGE

ISO 9001:2015  
9175.ING2