INGEZIM PESTIVIRUS COMPAC

Prod Ref: 12.BVD.K3

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I. TECHNICAL BASIS

FOR SERUM SAMPLES

The kit is based on a competitive enzymatic immunoassay (Blocking Elisa). We make a brief description of the technique below:

The antigen is fixed in a solid support (polystyrene plate). Each samples to be assayed will be added in one well. In case it has antibodies specific of pestivirus, they will bind to the antigen coating the plate and blocking the specific antigenic sites. If the sample does not contain specific antibodies, the specific antigenic site will stay free. After incubation with serum samples, we add a pestivirus p80/p125 specific monoclonal antibody (peroxidase conjugated). If the serum sample contains antibodies against pestivirus, they will not allow the conjugate bind to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate.

After washing the plate to eliminate all non fixed material, we can detect the presence or absence of labelled Mab by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

Hence, if there is a coloured well it means that the Mab conjugated has bound to the antigen because there are no antibodies blocking the antigen, so the sample could be marked as negative to Pestivirus. On the contrary, when the well does not give this colorimetric reaction means there are pestivirus specific antibodies blocking the antigen. In this case the sample must be marked as Positive to pestivirus.

FOR MILK & WHEY SAMPLES

The kit is based on an indirect enzymatic immunoassay (Indirect Elisa). We make a brief description of the technique below:

The antigen is fixed in a solid support (polystyrene plate). Each samples to be assayed will be added in one well. In case it has antibodies specific of pestivirus, they will bind to the antigen coating the plate. After incubation with milk samples, we add a bovine IgG specific monoclonal antibody (peroxidase conjugated). If the milk sample contains antibodies against pestivirus, the conjugate bind to them and will be fixed to the well.

After washing the plate to eliminate all not fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

Hence, if there is a coloured well it means that the conjugate has bound to the positive milk antibodies. On the contrary, when the well does not give this colorimetric reaction means there are not pestivirus specific antibodies.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to aliquot and store them at -20°C. Controls milk must be stored always at 4°C

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well..
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

SERUM SAMPLES

The sera samples need to be at 1/5 dilution in serum diluent to be tested (20 µl of serum + 80 µl of diluent for serum is the necessary quantity for one well: 100 µl in total)

If you wish to run titration, dilutions of the sera could be done (1/5, 1/10, 1/20, etc.).

MILK SAMPLES (individual and tank)

Must be tested undiluted. Fresh, refrigerated or previously frozen milk may be tested. To eliminate interference from lipids, partially skimmed milk should be either by centrifugation (15 min at 2000xg) or leaving overnight at 4°C to observe the formation of a lipid layer on the surface. Sample should be collected under the layer of lipids and that will be the one used in the assay. To obtain whey, freeze/thaw the milk sample.

VI. PREPARATION OF REAGENTS

• *Washing solution:*

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and 960 ml of water). When ready this solution remains stable between +2°C and +8°C.

• *Controls:*

FOR SERUM: Control sera must be diluted as sera samples

FOR MILK: Control milk must be tested as milk sample (undiluted).

• *Preparation of the conjugate (to make immediately before use):*

Dilute the needed quantity of conjugate provided in the Kit 1/100 with diluent

FOR SERUM use the conjugate **with yellow top**

- The necessary and sufficient quantity of conjugate for a complete plate is 60 µl of Conjugate in 6 ml of diluent.
- The necessary and sufficient quantity of Conjugate for an eight wells strip is 5 µl of conjugate in 0,5 ml of diluent.

FOR MILK use the conjugate **with green top**

- The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of Conjugate in 11 ml of diluent.
- The necessary and sufficient quantity of Conjugate for an eight wells strip is 10 µl of conjugate in 1 ml of diluent.

Shake very well the solution before use.

Prepare only the quantity needed for each time because the remainder volume has to be rejected

INSTRUCTIONS FOR SERUM SAMPLES TESTING

VII. TEST PROCEDURE

- All reagents (except conjugate) must be allowed to come to room temperature before use.
- Add 100 µl of Positive control serum, Negative control serum and samples (diluted as it's specified in the point V) to each well of the plate. We recommend running samples and control in duplicate. Incubate the plate **one-hour at 37 ± 1°C**.
- WITHOUT REMOVING SERUM SAMPLES**, add 50 µl of conjugate (prepared following previous instructions) to each well. Shake carefully the plate to allow the correct homogenisation of the components mixed on each well. Incubate the plate for **1 hour at room temperature (+25°C)**.
- Wash 5 times following the described procedure.
- Add 100 µl of substrate solution, to each well. Keep the plate for **15 min at room temperature**, in a dark place.
- Add 100 µl of stop solution to each well.
- Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution

VIII. .READING AND RESULT INTERPRETATION

VALIDATION OF THE TEST

The test could be considered valid if:

- The ODs for the Positive Control < 0.4
- The ODs for the Negative Control > 0.8

INTERPRETATION OF THE RESULTS:

Cut Off calculation:

- ⇒ *Positive Cut Off* = Negative control OD x 0,5
- ⇒ *Negative Cut Off* = Negative control OD X 0,55

Results Interpretation:

When you are running duplicate samples, OD values will be calculate for each sample as the arithmetic mean of both values.

- ⇒ All samples with OD higher than negative CUT OFF value must be considered as Negatives.
- ⇒ All samples with OD values lower than positive CUT OFF values will be considered as Positives
- ⇒ Samples with OD values between both cut-offs must be considered as doubt full. For these samples a new assay is recommended after 3 weeks.
- ⇒ In case of titration, the titter of the sample will be the higher dilution with an OD value lower than positive cut off.

EXAMPLE:

- OD of the Negative Control = 1,300
- OD of the Positive Control = 0,220
- Positive Cut Off = 1.300 x 0.5 = 0.650
- Negative Cut Off = 1.300 x 0.55 = 0.715

DILUTIONS	SAMPLE 1	SAMPLE 2	SAMPLE 3
1/5	0.163	0.291	1.235
1/10	0.176	0.412	1.284
1/20	0.211	0.567	1.401
1/40	0.325	0.842	1.180
1/80	0.505	1.055	1.320
1/460	0.617	1.330	1.324
1/320	0.897	1.284	1.408
1/640	1.222	1.357	1.260

In this example samples 1 and 2 are positive with titters 1/160 and 1/20 respectively. Sample 3 is negative.

INSTRUCTIONS FOR MILK AND WHEY SAMPLES TESTING

VII. TEST PROCEDURE

1. Add 100 µl of samples milk and Positive and negative control milk (as it's specified in the point V and VI) to each well of the plate. We recommend running samples and control in duplicate. Incubate the plate **45 min at 37 ± 1°C**.
2. Wash 4 times following the described procedure
3. Add 100 µl of conjugate (prepared following previous instructions) to each well. Incubate the plate for **30 min at 37 ± 1°C**.
4. Wash 5 times following the described procedure.
5. Add 100 µl of substrate solution, to each well. Keep the plate for **10 min at room temperature**, in a dark place.
6. Add 100 µl of stop solution to each well.
7. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution

VIII. .READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

VALIDATION OF THE TEST

The test could be considered valid if:

- Abs (Control +) > 0,8
- Abs (Control -) < 0,24

INTERPRETATION OF THE RESULTS:

When you are running duplicate samples, OD values will be calculate for each sample as the arithmetic mean of both values.

S/P calculation:

$$S/P = OD \text{ sample} / OD \text{ control+}$$

Results Interpretation:

⇒ $S/P > 0.3$ sample is positive

⇒ $S/P \leq 0.3$ sample is negative

INGEZIM PESTIVIRUS
12.BVD.K3

ENGLISH

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: