

## INGEZIM PESTIVIRUS Compac

R.12.BVD.K3

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

## INgezim PESTIVIRUS Compac

R.12.BVD.K3

**INgezim PESTIVIRUS Compac** is an enzymatic assay based on a competition ELISA technique, which uses a monoclonal antibody (MAb) specific for ruminant pestivirus (BVDV & BDV) p80 protein, and a recombinant antigen (pestivirus p80 protein).

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with BVD antigen (recombinant protein p80). Serum samples are added and incubated.
2. If the samples contain specific antibodies to the p80 protein of pestiviruses, they will bind to the antigen.
3. When a MAb-PO specific of p80 protein is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected animals or vaccinated animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection of specific antibodies to the highly conserved non-structural NS3 protein (p80) of pestiviruses, in bovine, ovine, goat serum and milk samples.

### INTERPRETATION OF THE RESULTS

Two cut offs are used for the results interpretation: positive and negative. The samples will be considered as **Positive** when their OD value is equal to or lower than the positive cut off. The samples will be considered as **Negative** when their OD value is equal to or higher than the negative cut off. The samples with an OD value between both cut offs will be considered as **Doubtful**.

## VALIDATION

### Study 1: Correspondence with SN (seroneutralization)

#### Cattle:

**Infected animals:** 73 bovine sera belonging to 4 herds infected with BVDV, free of PIs and not vaccinated were analyzed. Results obtained indicate that the correspondence between INgezim® PESTIVIRUS COMPAC and SN was **92.5% (90% sensitivity & 96% specificity)**.

**Vaccinated animals:** These 73 cattle were vaccinated with an inactivated vaccine (Bovilis® BVD: containing mainly gp53) at days 0 (1<sup>st</sup> dose) & 30 (2<sup>nd</sup> dose). At days 0, 30 & 60 animals were bled. Results obtained showed **82%** correspondence between both techniques.

- **Small ruminants:** 57 sheep and 48 goats belonging to 8 and 9 sheep and goats flocks were analysed by SN & INgezim Pestivirus Compac. Results obtained showed a correspondence between both techniques higher than 99.9%.


### Study 2: Pestivirus seroprevalence studies in small ruminants.

In total, 201 sheep, 112 goat and 147 mixed flocks were sampled. Eighteen per cent (497 out of 2,797), 6% (85 out of 1,520) and 21% (763 out of 3,609) of the sera from sheep, goat and mixed flocks, respectively, were positive. This study indicates that INGEZIM PESTIVIRUS Compac showed the best balance in sensitivity-specificity, as well as the lowest amount of undetermined results for the detection of antibodies to pestiviruses in sheep sera compared with other three commercially available ELISAs using the virus neutralisation test (VNT) with the cytopathic (cp) BVDV Oregon C24V strain as reference test (Valdazo et Álvarez, 2003. 6th International Congress of Veterinary Virology, Saint Malo, Francia p.134)

### Study 3: Use of milk samples.

524 tanks samples previously classified by Svanovir® BVDV-Ab were analyzed. Correspondence between both assays was 88%.

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### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB)

REGISTRATION NUMBER 1054 RD  
PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**  
Stored at 2°C-8°C

Ed.020218