INGEZIM AGID

Prod Ref: 12.BB2.K7

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

I. TECHNICAL BASIS

The AGID technique allows to detect presence of antibodies specific of Brucella in serum. This technique is based on the use of an antigen containing mainly lipopolysaccharide (PLS) and another polysaccharide named HN (both of them present in the cellular membrane). This antigen and the antibodies present in a positive serum diffuse through the

agarosa gel to form a complex. When the complex is formed it precipitates and the reaction is manifested by the development of a visible line of precipitate in the interface. Due to the use of these 2 antigens, the assay allows to distinguish between infected and vaccinated animals.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Do not mix reagents or use instructions from different kits.
3. Avoid any contamination of the kit reagents.
4. Do not use components after expiry dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Systematically include a positive control and a negative serum each time the assay is run.

III. STORAGE OF COMPONENTS

The agar from the plate may be altered depending on the storage conditions. It is therefore recommended to store them in conditions that prevent desiccation: always keep the plates well closed and upside down (to prevent condensation drops in the lid from falling onto the agar layer), and at

temperatures between +2 ° C and +8 ° C.

The positive and negative controls and antigen should be kept frozen. For partial uses of the kit, it is convenient to aliquot these reagents in order to avoid their continuous freezing and thawing.

IV. TEST PROCEDURE


1. Dispense 15 µl of positive control and 15 µl of negative control into 2 lateral wells.
2. Dispense 15 µl of each serum in the rest lateral wells
3. Dispense 15 µl of antigen into the central well.
4. The thus inoculated plates are incubated in a humid chamber at room temperature. Incubate the plate at room temperature for 24-48 hours

V. LECTURE AND RESULT INTERRETATION

The results are read after 24 and 48 hours, with oblique illumination.

A. Interpretation of the results:

- The absence of precipitation lines means that the animal is not infected.
- The presence of precipitation lines can mean:
 - a. **Double band** (the band corresponding to the NH antigen is observed closer to the serum well, and the band corresponding to the LPS antigen, near the central well). It means that the animal is infected and is probably excreting the bacteria. Some very recently vaccinated animals (two months or less after vaccination) could also produce that double band of precipitation, characteristic of infected animals.
 - b. **A single band** corresponding to the LPS antigen. This reaction occurs transiently in vaccinated animals, although it could occur in a relatively low proportion of infected animals (which react only with the LPS band instead of with both antigens, LPS and NH).

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

IT-73840
IT-73780

ISO 14001:2015
9191.INGE

ISO 9001:2015
9175.ING2

Distributed in

by: