## INgezim Total Milk CROM

Prod. Ref: 30.MLK.K42

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## INgezim Total Milk CROM

Prod. Ref: 30.MLK.K42

Direct immunochromatographic assay, for specific detection of milk protein in food and environmental samples

### KIT COMPONENTS:

| Component  | Units |    |
|--|-------|----|
|  |       |    |
| INgezim Total Milk CROM assay strips   | 10    | 25 |
| Bottle containing 60 ml of extraction solution, ready to use                         | 1     | 3  |
| Plastic disposable pipettes  | 10    | 25 |
| Empty tubes for extraction procedure   | 10    | 25 |
| Microtiter 8-well strips for dispensing the extracted samples and running the assay. | 2     | 4  |

### OTHER NEEDED MATERIALS NOT PROVIDED WITH THE KIT:

- Pipette for measuring 5ml
- Grinder, stomacher, mortar or any other, manual or automatic homogenization system for an adequate sample crushing.

### STORAGE CONDITIONS OF THE KIT:

All components should be stored between +4 °C and +25 °C, always inside its original packaging till the time of use.

### TEST BASIS:

The test is based on a lateral flow assay, a migration technique which uses specific antibodies to casein and  $\beta$ -lactoglobulin, and it is able to detect residues of these proteins in all food matrix and environmental samples (there is a specific protocol for this kind of samples. Please consult the manufacturer).

### DETECTION LIMITS:

- INgezim Total Milk CROM is a qualitative assay. In case of interest in casein and/or  $\beta$ -lactoglobulin quantification (ppm), we recommend to use the ELISA INgezim Caseina assay (Prod. Ref. 30.CAS.K2) and/or ELISA INgezim  $\beta$ -lactoglobulina (REF 30.BLG.K2).
- The detection limit of the INgezim Total Milk CROM is 0.25 ppm of milk protein diluted in the extraction buffer.
- INgezim Total Milk CROM is able to detect up to 0.25 ppm of milk protein in soy milk matrix fortified with milk protein. Considering that the sample to be assayed with the kit has been diluted 1/10 in extraction buffer, sample should contain more than 2.5 ppm of casein to obtain a positive result.
- Samples very viscous, dense or with a high content of fat could migrate incorrectly along the membrane. Thus, signal can be affected (attenuation or disappearance of test and control line).
- The test INgezim Total Milk CROM is specific for casein and  $\beta$ -lactoglobulin, and do not cross react or recognises other milk proteins.
- When it is necessary to differentiate if the milk contamination comes from casein or from  $\beta$ -lactoglobulin, we recommend to use the INgezim Caseina CROM assay (30.CAS.K42) and/or INgezim  $\beta$ -lactoglobulina CROM assay (30.BLG.K42).

### WARNINGS

1. Read the instructions carefully.
2. Bring all reagents to room temperature before use.
3. Do not mix reagents from different kits or batches.
4. Avoid contamination of reagents.
5. Do not use kits after the expiration date.
6. Do not eat, drink or smoke while handling samples and / or reagents.

(\*)NOTE: The more sample amount is used for the extraction, the more representative will be the analysis and the more reliable the result. If you want to extract more amount of sample, be sure to always maintain the relationship **1:10** (sample: extraction buffer).

INMUNOLOGIA Y GENETICA APLICADA, SA (INGENASA)

C/ Hnos Garcia Noblejas, 39 – 28037 MADRID (Spain)

Tel: (+34)91 3680501

www.ingenasa.com

Version doc:150218

### INSTRUCTIONS:

#### A.- SAMPLE PREPARATION (\*):

##### Solid Samples

1. Crush the sample to obtain a powder consistency as thin as possible. Use a mortar or a mechanical homogenizer. If you are going to analyze more than one sample, wash properly the homogenizer with soap, rinse thoroughly with water and finally dry it with 70 % ethanol.

##### Liquid Samples

1. Shake the sample to ensure it is homogeneous and that you are taking a representative part of it.
2. Weigh **0.5 g** of sample and place it into a provided tube. Add **5 ml** of supplied extraction buffer.
2. Take **0.5 ml** of sample and place it into a provided tube. Add **4.5 ml** of supplied extraction buffer.
3. Stir for 15-30 seconds using a vortex mixer or agitator to ensure homogenization.
4. In case of solid samples, allow the sample to settle for a few minutes.
5. Add 10 drops of the supernatant in a clean well (8-well strips provided). For samples with high fat content, avoid taking the fat layer of the supernatant.


#### B.- TEST PROCEDURE:

Open the tube containing the strips, JUST before performing the assay, and close it immediately.

Take the strips and put them vertically, one in each well containing the prepared samples. The white end of the strip should be pointing at the sample direction.

Keep there during 10 minutes. After these 10 minutes, take the strips out and put them on a flat and (preferably) white surface, for reading the results.

#### C.- READING AND INTERPRETATION OF THE RESULTS:

**POSITIVE NEGATIVE NOT VALID**

• **POSITIVE:** Two red lines appear on the reading zone, being the upper one the control band and the bottom one the results band. The color intensity of the results line may vary, but it could be not proportional to the milk protein concentration in the sample.

• **NEGATIVE:** Only the upper red line appears on the reading zone.

• **NOT VALID:** There are two possible situations: 1) no lines are visible or 2) only the bottom red line appears. Not valid results may be caused by damaged reagents or by an incorrect manipulation. In this case, the test should be repeated using a new strip.

*We recommend not to consider any sample as negative until 10 minutes have passed.*