
## INgezim® TOTAL MILK CROM

R.30.MLK.K42

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
## INgezim® TOTAL MILK CROM

R.30.MLK.K42

**INgezim® TOTAL MILK CROM** is a direct immunochromatographic assay, based on Lateral Flow technique, which uses 2 monoclonal antibodies (MAbs) specific of two different epitopes of Casein (CAS) and a specific  $\beta$ -lactoglobulin ( $\beta$ -Lac) polyclonal antibody.

### TECHNICAL BASIS OF THE KIT

The device consists on a strip of nitrocellulose (NC) coated with:

- **Lower end:** Contains a mix of MAb specific of casein and a polyclonal specific of  $\beta$ -Lac conjugated to latex particles.
- **Test Line:** Contains a mix of MAb specific of casein and a polyclonal specific of  $\beta$ -Lac.
- **Control line:** Contains a MAb specific to the control protein

The strips are introduced in a well containing the prepared samples.

If there are either of these proteins in the sample, they will bind to the Ac conjugated to latex particles. The complex Ag-Ac-Latex will flow through the membrane until be captured by the Ab adsorbed in the Test Line.

Control protein bound to latex will be captured in the Control Line.

### APPLICATION

Detection of residues of Casein and  $\beta$ -Lactoglobuline in all food matrix and environmental samples. The kit does not differentiate between both kind of proteins.

### INTERPRETATION OF THE RESULTS

If the sample is positive to CAS or to  $\beta$ -Lac, two red lines will appear on the reading zone, being the upper one the control band and the bottom one the results band. If the sample is negative to both proteins, only the upper red line will appears on the reading zone.

### VALIDACIÓN

#### Analytical sensitivity

- Able to detect up to 0.25 ppm of milk protein in soy milk matrix fortified with milk protein. Considering that the sample to be assayed with the kit has been diluted 1/10 in extraction buffer, sample should contain more than 2.5 ppm of casein to obtain a positive result.

#### Diagnostic sensitivity and specificity

91 samples (liquid and solid) have been analyzed. The evaluated matrixes consist in: processed foods based on cereals, alcoholic beverages (liquors, wine), non-alcoholic beverages (juices, soft drinks, smoothies, etc.), baby food, freeze-dried foods, bakery products, cookies, industrial snacks, sauces for seasoning, meat preparations and sausages, chocolate, legumes (it has been described cross-reaction in some references and with commercial tests). Results were compared with the obtained using the specific tests to each protein.

		INGEZIM CASEIN / $\beta$ -Lac ROM	
		POS	NEG
TOTAL MILK CROM	POS	43	5*
	NEG	0	43

\*3 of the 5 samples refer milk in the label.

**Correlation with the specific assays is 94.5%**

### COMPOSITION OF THE KIT

Ingezim Total Milk CROM assay strips
Bottle containing 60 ml of extraction solution
Plastic disposable pipettes
Empty tubes for extraction procedure
Microtiter 8-well strips.

PRODUCT MANUFACTURED BY INGENASA

**SHELF LIFE: 24 months**  
Stored: 4°C-25°C

Ed.080818