INGEZIM PPA COMPAC

Prod Ref: 11.PPA.K3

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I. TECHNICAL BASIS

This kit is based on a blocking enzymatic immunoassay (Blocking Elisa). We make a brief description of the technique below:

The antigen is fixed in a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen absorbed on plate while if the serum sample does not contain specific antibodies they will not bind the antigen. If we add a specific monoclonal antibody (Mab) against the viral antigen coated to the plate (conjugated with peroxidase), it will compete with the antibodies of the serum. If the serum samples contains specific antibodies, they will

not allow the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material from the plate, we can detect the presence or absence of labelled Mab by adding the specific substrate that in presence of the peroxidase will develop a colorimetric reaction.

The antigen coated to the plate in our kit consists of a purified protein extract from the virus (the VP72), which is the major structural protein from the ASFV and the most antigenic one.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
10. Substrate must be handle with care, it is very sensible to light and contamination.
11. Stop solution is a strong acid. Handle with care.

III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C.

Once opened, control sera are stable for one month between +2°C and +8°C. In case that they are not going to be used in this period, we recommend to store them under -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.

- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
- Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLE

Samples to be used in the assay must be porcine sera. **Do not use highly haemolysed or contaminated samples. You might obtain false positive results.**

To be assayed a ½ dilution must be done in the supplied serum diluent (DE31).

This dilution could be done directly on the plate by adding 50µl of serum diluent and 50 µl of sera into each well.

VI. PREPARATION OF REAGENTS

• **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (i.e. 40 ml of concentrate plus 960 ml of dH2O). When ready this solution remains stable between +2°C and +8°C until the expire date described at the label of the concentrated solution.

• **Control sera:**

Once opened, control sera are stable for one month between +2°C and +8°C. In case that they are not going to be used in this period, we recommend to store them under-20°C.

The control sera need to be diluted ½ in diluent prior to be used, in the same way than the sera samples (you can do the dilution directly into the plate by adding 50 µl of diluent plus 50 µl of control sera).

• **Preparation of the conjugate: to make immediately before use.**

Dilute the needed quantity of conjugate provided in the kit 1/100 into the supplied conjugate diluent (DE01) (e.g. 110µl of conjugate into 11 ml of conjugate diluent is enough quantity for one plate).

Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature (20-25°C) before use.
2. Add 50 µl of supplied serum diluent to each well. Add 50 µl of positive control sera to two wells (e.g. A1 and B1), and 50 µl of the negative control sera (e.g. A2 and B2). Add 50 µl of sera samples to test on each remainder wells. We recommend the use of two wells for each control. Seal the plate and **incubate for 1 h±5min at 36±1°C or overnight (16-20 hours) at 20-25°C.**
3. Empty the wells into a container with NaOH solution and wash 4 times as specified in point IV.
4. Add 100 µl of specific conjugate (prepared following previous instructions) to each well. Seal the plate and **incubate for 30 minutes at 36±1°C.**
5. Wash 5 times as previously described
6. Add 100 µl of substrate to each well. Keep the plate for 15 min at room temperature (20-25°C).
7. Add 100 µl of stop solution, to each well.
8. Read the OD of each well at 450nm within 5 min after the addition of stop solution

VIII. READING AND RESULT INTERPRETATION

Validation of the test

The test could be considered valid when the OD of the Negative Control (NC) is, at least ,4 times higher than the OD of the Positive Control (PC):

$$NC / PC \geq 4$$

Interpretation:

A. - Cut off calculation:

$$\text{Positive Cut Off} = NC - [(NC - PC) \times 0.5]$$

$$\text{Negative Cut Off} = NC - [(NC - PC) \times 0.4]$$

Where: NC = OD of Negative Control Serum
PC = OD of Positive Control Sera.

For calculating the blocking % (x %) of a sample:

$$X\% = \frac{NC - \text{SAMPLE OD}}{NC - PC} \times 100$$

≥ 50% positive

≤ 40% negative

Between both values, doubtful

B. - Results Interpretation:

When you are running duplicates, OD of the sample, will be calculated as the arithmetic mean of OD values in both wells.

- Serum samples with an OD lower than Positive Cut Off, are considered as POSITIVES to ASFV antibodies.
- Serum samples with an OD higher than Negative Cut Off, will be considered NEGATIVE to ASFV antibodies.
- Serum samples with OD values between both cut off are considered as doubtful. We recommend re-testing these animals one more time or applying a different technique to check this serum (IFI, Indirect ELISA, etc.).

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: