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
I. TECHNICAL BASIS

The reagents will be used a Double antibody sandwich enzyme immunoassay.

Here we describe the technical basis of the assay: a specific monoclonal antibody (Mab) against Rabies virus glycoprotein is coated on a solid support (polystyrene plate). When a sample containing viral antigen is added, the Mab will capture antigen. After washing (to remove all material not bound to the plate), we add a second specific MAb to glycoprotein labelled with biotin (conjugated I).

This MAb will bind to it, and after incubation with a second conjugate (streptavidin-peroxidase), we could demonstrate the presence of this conjugate (conjugate II) by adding a chromogenic substrate that in presence of the enzyme will produce a colorimetric reaction.

II. PRECAUTIONS:

1. Read the instructions of use carefully.
2. **Bring all reagents to room temperature (20°-25°C) prior to use.**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. Stop solution is a strong acid. Handle with care.
9. Substrate is very sensitive to light and contamination.

III. STORAGE:

Store all reagents between +2°C y +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

1. Throw out the content of out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
2. Dispense a volume of 300 µl of washing solution on each well.
3. Shake delicately the plate, avoiding the contamination between wells.
4. Turn over the plate brusquely to empty the wells.
5. Repeat the process as much times as is indicated on the instructions of the kit.
6. Prior to empty the content of the last washing step, verify that the next reagent to be added to the After the last step of washing shake the plate turned over an absorbent filter paper

V. PREPARATION OF SAMPLES:

To titrate the samples it is recommended to make eight dilutions in fold of 2 from 1/10 dilution in sample diluent I.

VI. PREPARATION OF REAGENTS

◆ **Washing solution:**

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C.

◆ **Preparation of the conjugate I: to make immediately before use.**

Dilute 1/100 the needed quantity of conjugate with diluent I (e.g. 10 µl of conjugate I with 1 ml of diluent I is more than enough quantity for a strip of 8 wells). Shake very well the

solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

◆ **Preparation of the conjugate II: to make immediately before use.**

Dilute 1/100 the needed quantity of conjugate II into diluent (e.g. 10 µl of conjugate II with 1 ml of diluent II, is enough quantity for a strip of 8 wells). Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected

VII. TEST PROCEDURE:

1. Take out the plate (or strips to be used) and maintain 15 min. to room temperature before starting the test.
2. Add 100 µl/ pocillo of each sample. Seal the plate and **incubate for 1 h a 37°C.**
3. Throw away the liquid contained in the plate on a recipient containing NaOH 0.1 M. Wash 3 times following the procedure previously described.
4. Add 100 µl of conjugate I, prepared following the previous instructions, to each well. Seal the plate and **incubate for 1 h a t 37°C.**
5. Wash 4 times following the described procedure.
6. Add 100 µl of Conjugate II, prepared following the previous instruction, to each well. **Incubate 30 minutes at room temperature (25°C).**
7. Wash 5 times more following the described procedure.
8. Add 100 µl of substrate. **Keep the plate for 10 min at room temperature.**
9. Add 100 µl of stop solution to each well.
10. Read the absorbance of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. INTERPRETATION OF THE RESULTS

1. In this product, controls are not included, thus the users must prepare their own standards and establish validation and interpretation criteria in agreement with their process.
2. The product has not been evaluated to be used for diagnosis.
3. INGENASA is not responsible of the results derived from an inadequate use of the reagents or for the wrong selection of standards.

ATTENTION! This product is not a validated one. Just for internal uses. The use of EIA could be prohibited in some countries. Please, observe local state and federal rules for proper use.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by:

IT-73840
IT-73780

ISO 14001:2015
9191.INGE

ISO 9001:2015
9175.ING2