INgezim TB CROM Ab

Prod Ref: 11.TBP.K41

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NEGATIVO

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I. TECHNICAL BASIS

Ingezim® TB CROM Ab is a lateral flow assay that uses a protein of *Mycobacterium bovis* to detect specific antibodies in serum and blood samples.

The test is based on the use of two coloured latex microspheres: red microspheres covalently bound to MPB83 protein of *M. bovis*, and green microspheres bound to a control protein that are used as test control.

The membrane presents a test line (T), in which the MPB83 protein of *M. bovis* was dispensed; and a control line (C), containing a specific MAb for the control protein.

If the specific antibodies are contained in the sample, they will react with the latex particles which are coated with MPB83. The red *latex particles/protein/antibodies* complex will flow through the membrane and bind to the protein contained in the test line, where a colour line will appear. If specific antibodies are not present in the sample, no colour line will appear.

II. STORAGE

All the reagents must be kept between +4°C and +25 °C.

III. PRECAUTIONS Read the use instructions carefully.

1. Bring all reagents to room temperature (20°-25°C) prior to use.
2. Do not mix instructions or reagents from different kits.
3. Avoid any contamination of the reagents of the kit.
4. Do not use components after expiration date and do not mix components from different batches.
5. Do not eat, drink, or smoke where specimens or kit reagents are being handle.

IV. SAMPLES


The test has been designed for serum and blood samples. Serum samples must be fresh. Blood samples must be fresh or refrigerated up to 4 days at 2-8 °C, and collected with any anticoagulant (EDTA, heparin, etc.). Blood samples without any anticoagulant could contain micro-clots that can block the device and result in nonspecific reactions.

It is very important to bring samples to room temperature (20-25 °C) prior to use.

V. TEST PROCEDURE

- **Allow all kit components and samples to reach room temperature, before use.**
- Remove the device from the aluminium pouch, place it on a flat surface and discard the desiccant.
- Take **the serum/blood sample** using the pipette and add it to the bottle that contains running buffer. Insert the dropper in the bottle.
- Shake the dropper bottle to assure a homogeneous solution is obtained. Add **4 drops** of the solution onto the sample window (it is very important to add the sample dilution drop by drop in order to have a good lateral flow speed).
- Wait **12 minutes** and read the results.

IT IS VERY IMPORTANT TO FOLLOW THESE INSTRUCTIONS AND RESPECT THE INDICATED VOLUMES TO OBTAIN VALID RESULTS

1. Add the sample to the bottle, insert dropper and shake
2. Add 4 drops to the strip and read after 12 min.


VI. INTERPRETATION OF RESULTS

Interpretation of results should be done as follows, 12 minutes after the addition of the sample:

POSITIVE result: Two lines appear on the reading zone: a red line (positive) and a green line (control). The colour intensity of the red line may vary depending on the antibody concentration.

NEGATIVE result: Only one green line appears on the reading zone.

INVALID result: There are two possible situations: 1) no lines are shown on the reading zone; or 2) only one red line appears. Not valid results may be caused by damaged reagents or by an incorrect manipulation. In this case, the test should be repeated using a new strip.

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Distributed in _____ by:

Empty box for distributor information.

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INMUNOLOGIA Y GENETICA APLICADA, S.A.