INgezim FCoV IC

Prod. Ref: **1.6.FCV.K41**

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INgezim FCoV IC

Prod. Ref: **1.6.FCV.K41**

Immunocromatographic test for feline coronavirus (FCoV) specific antibodies detection in cat serum, blood or plasma samples.

INMUNOLOGIA Y GENETICA APLICADA, SA (INGENASA)
C/ Hnos Garcia Noblejas, 39 – 28037 MADRID (Spain)

Tel: (+34)91 3680501
www.ingenasa.com

Versión doc: 022016

KITS COMPONENTS:

Component	Units
Immunocromatographic devices	12
Dropper containing chromatographic buffer	1
Micropipette for 10µl	1
Micropipette teps	12

OTHER MATERIALS OR REAGENSTS NEEDED NOT PROVIDED:

– None

KIT CONSERVATION:

All components of the kit should be stored at + 4 ° C and + 25 ° C, in its original packaging until use.

PRINCIPLE OF THE TEST:

The test is based on the indirect immunocromatographic technique, and has been designed to detect antibodies to FCoV in cat serum, blood or plasma samples.

PRECAUTIONS:

- Carefully follow the instructions.
- Place all reagents to room temperature at least 20 min before use.
- Never open the foil pouch containing the device if not used immediately. It could be hydrated and become unusable. Open only immediately before use.
- Do not mix instructions or reagents from different kits.
- Avoid any contamination of reagents.
- Do not use the kit once exceeded their expiration date.
- Do not eat, drink or smoke when samples and / or reagents are handled.
- Do not touch the chromatographic membrane with any object, especially with the dropper tip.

INSTRUCTIONS:

A.- SAMPLE PREPARATION:

The sample of choice for diagnostic purposes, can be serum, plasma or whole blood with anticoagulant (EDTA, heparin ...).

The test should be performed after extraction of the sample. If this is not possible, samples should be stored at 2-7 ° C for up to 5 days. If test is planned to be done after more than 5 days, samples should be stored frozen.

B.- TEST PROCEDURE:

Remove from the bag the diagnostic device and place it on a flat surface. NOTE: Only open the bag at the time of going to test.

Sample addition to the device: Using a micropipette, add in the "S" window 10 µl sample (serum, plasma or whole blood).

CAUTION: The pipette must be placed vertically on the window and keeping the tip 1 cm to avoid touching the membrane.


Add 3 drops of chromatographic buffer in the "S" window.

C.- READING AND RESULTS INTERPRETATION:

Keep the device for a flat surface and untouched for 10 minutes, after which proceed to reading the result:

POSITIVE
result

NEGATIVE
result

It is advisable not to give any negative sample definitely until passed 10 minutes after placing the sample on the device.