
## INgezim HAEMOPHILUS PARASUIS

R.11.HPS.K1

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
## INgezim HAEMOPHILUS PARASUIS

R.11.HPS.K1

**INgezim Haemophilus parasuis** is based on the Indirect ELISA technique which uses a polyclonal antibody specific of porcine immunoglobulines as conjugate and a crude extract of proteins of Haemophilus parasuis (HPS) as antigen.

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with a crude extract of HPS. Serum samples are added and incubated.
2. If the samples contain specific antibodies to Haemophilus parasuis, they will bind to the antigen.
3. When the conjugate specific of porcine Igs is added, it will bind to the Igs bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

INgezim Haemophilus parasuis has been designed for diagnosis of infections by HPS (Glasser Disease) detecting increase of specific antibodies (seroconversion) and for evaluation of the vaccination response

### INTERPRETATION OF THE RESULTS

The assay uses two Cut offs: Samples showing OD values higher than or equal to the positive cut off will be considered **positive**. Samples showing OD values lower than the negative Cut off must be considered **negative**. There is a range of doubtful samples between these two values.


### VALIDATION

#### 1. Early detection of specific antibodies

The experiment was conducted with pigs originating from a high health status herd (PCR, isolation). 30 pigs (group A) were transferred to an isolated barn when about 5 month old. They were bled on days 1, 8 and 17 and the sera were examined using INgezim HPS. On day 24 pigs from another herd HPS positive were introduced and mixed with the pigs already present. They were bled on day 38 and sera were examined using INgezim HPS. These pigs demonstrated low (<0.2) and stable S/P from day 1 to 24. No clinical signs were observed and bacterial isolation and PCR performed on day 8 and 17 were negative for HPS suggesting that they were not infected. A few days after mixing, pigs from group A demonstrated clinical signs compatible with Glasser's disease (fever, anorexia, lameness). Some of them died and demonstrated lesions compatible with Glasser's disease. HPS was isolated from several pigs (untypable strain). Surviving pigs demonstrated a significant increase in S/P on day 38 (**14 days after mixing with pigs from the HPS positive source**)

#### 2. Specificity

Determination of the specificity is very difficult due to the fact that HPS is present in most of herds including those with a high health status. There were used 50 animals from a HPS free herd, analysed in 4 separate occasions 3 weeks apart. All these animals were negative during the process.

— positive cut off — negative cut off

#### 3. Sensitivity to different serovars

Sera of animals experimentally infected with **serovars: 1, 4, 5, 7, 8, 9, 10, 11, 13, 14 and 15** were analysed. This assay is able to detect specific antibodies to all of them.

### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive control
- Vials with Negative control
- Vials with Conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA  
Registration number 1215 RD

IT-73840 ISO 9001:2015  
IT-73780 9191.INGE 9175.ING2

SHELF LIFE: **18 months**  
Store at 2°C-8°C

Ed.020217