hMPV

## Anticuerpo Monoclonal específico de MPVh (F) hMPV (F) specific Monoclonal Antibody

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*The human Metapneumovirus belongs to the Paramyxoviridae family, Pneumovirinae subfamily, to which also belongs VRS. There exist two different genetic lineages, each one consisting of two different subtypes. MPVh is an enveloped negative simple stranded RNA virus. The main proteins are: N (nucleocapsid), M (matrix), and 3 glycoproteins in the envelope (G, F, SH). This virus exclusively affects humans. All of the evidences obtained since it was discovered, suggest that it is one of the main etiologic agents causing low respiratory infection in children.*

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
*The hybridoma which produces the monoclonal antibody has been obtained by the fusion of lymphocytes from Balb/c mice's spleen with myeloma X63/Ag8653 cells. The IgG has been purified by affinity chromatography, showing a purity of 99%.*

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*Detection of F protein of human Metapneumovirus by immunofluorescence (IFA)*

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*Infected cells are deposited onto microscope slides, dried and fixed with acetone. The slides are incubated with the primary antibody (See recommended use conditions for dilution. Note that the supernatant is supplied ready to use) for 20 min in a wet chamber at 37°C, wash 3 times with phosphate-buffered saline (PBS), and then incubate with an anti-mouse immunoglobulin conjugated with fluorescein-isothiocyanate at a dilution of 1/100 in PBS with Evans Blue (0.01%) for 20 min at 37°C. The slides are washed and dried and analyzed using a fluorescence microscope.*

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